

Optical Systems

Table of Contents




F-Theta Lens
page 5.3



Compact Beam Expander
page 5.4



Zoom Beam Expander
page 5.4



Top Hat Beam Shaping Lens – FBS
page 5.5



Continuously Variable Attenuator / Beamsplitter 990-0060
page 5.8



Variable Attenuators for linearly polarized laser beam 990-0070
page 5.9



Motorized Variable Attenuator for linearly polarized laser beam 990-0070M
page 5.10



Variable Attenuators for linearly polarized laser beam 990-0071
page 5.12



Motorized Variable Attenuator for linearly polarized laser beam 990-0071M
page 5.13



Variable Attenuator for femtosecond laser pulses 990-0072
page 5.14



Variable Attenuator for femtosecond and Nd:YAG laser pulses 990-0073
page 5.15



Compact Motorized Laser Power Attenuator 990-0075
page 5.16



Compact Variable Laser Power Attenuator **990-0076**
page 5.17



Precision Spatial Filter **990-1000**
Find more at EksmaOptics.com



Y-Z Positioner for lens, pinholes and objectives **990-0100, 990-0200**
Find more at EksmaOptics.com



Y-Z Positioners for lens, pinholes and objectives **990-0050, 990-0051**
Find more at EksmaOptics.com



Precision Pinholes
Find more at EksmaOptics.com



Unmounted Iris Diaphragms
Find more at EksmaOptics.com



Mounted Iris Diaphragms
Find more at EksmaOptics.com



Mounts for iris diaphragms
Find more at EksmaOptics.com



Motorized Iris Diaphragms **995 Series**
Find more at EksmaOptics.com



Motorized Iris Diaphragms **996 Series**
Find more at EksmaOptics.com



Motorized Iris Diaphragms **997 Series**
Find more at EksmaOptics.com



Variable Wheel Attenuator **990-0604**
Find more at EksmaOptics.com



Closed Variable Wheel Attenuator **990-0704**
Find more at EksmaOptics.com



Filters Holder with 90° Flip **990-0400**
page 5.18



Motorized Variable Two Wheels Attenuators **991-0602**
Find more at EksmaOptics.com



Motorized Closed Variable Two Wheels Attenuators **991-0702**
Find more at EksmaOptics.com



Air-cooled Beam Dump **990-0800**
page 5.19




Water-cooled Beam Dump **990-0820**
page 5.19


F-THETA LENS




F-Theta lenses are designed to provide a flat field on the image plane for scanning and engraving applications where a high power laser and a set of rotating mirrors are used to scan across a given field.



A. Lens Diameter – 90 mm



B. Lens Diameter – 104 mm



C. Lens Diameter – 104 mm

BEST MIRROR PLACES m1/m2 – 16/16 mm, screw size – M85×1

Wavelength – 1064 nm, Lens Diameter – 90 mm

Focus length, mm	Working distance S, mm	Max. scan area, mm ²	Max. scan angle, θ max	Input beam diameter, mm	Spot size, μ m	Drawing	Catalogue number	Price, EUR
100	114	70×70	±28°	12	16	A	150-1001	420
160	232	140×140	±28°	12	26	A	150-1601	420
210	287	145×145	±28°	12	34	A	150-2101	420
254	284	175×175	±28°	12	31	A	150-2541	420
290	355	200×200	±28°	12	31	A	150-2901	420
330	465	220×220	±28°	12	40	A	150-3301	420
420	467	300×300	±28°	12	50	A	150-4201	420

Wavelength – 532 nm, Lens Diameter – 90 mm

Focus length, mm	Working distance S, mm	Max. scan area, mm ²	Max. scan angle, θ max	Input beam diameter, mm	Spot size, μ m	Drawing	Catalogue number	Price, EUR
100	114	70×70	±28°	12	16	A	150-1002	460
160	180	110×110	±28°	12	16	A	150-1602	460

Wavelength – 355 nm

Focus length, mm	Working distance S, mm	Max. scan area, mm ²	Max. scan angle, θ max	Input beam diameter, mm	Spot size, μ m	Drawing	Catalogue number	Price, EUR
100	136	70×70	±25°	7	10	B	150-1003D1	930
160	186	110×110	±25°	7	15	B	150-1603	930

BEST MIRROR PLACES m1/m2 – 24/24 mm, screw size – M85×1


Wavelength – 1064 nm, Lens Diameter – 104 mm

Focus length, mm	Working distance S, mm	Max. scan area, mm ²	Max. scan angle, θ max	Input beam diameter, mm	Spot size, μ m	Drawing	Catalogue number	Price, EUR
163	185	110×110	±28°	20	17	C	151-1631	520
210	255	150×150	±28°	20	24	C	151-2101	520
254	285	175×175	±28°	20	31	C	151-2541	520
420	467	300×300	±28°	20	55	C	151-4201	520
650	697	400×400	±25°	20	85	C	151-6501	520

COMPACT BEAM EXPANDER



A laser beam expander is designed to increase the diameter of a collimated input beam to a larger collimated output beam. EKSMA OPTICS offers compact Galilean type beam expanders for 1064 nm, 532 nm and 355 nm wavelengths. Compact beam expander has the possibility to be adjusted for the input beam divergence angle to obtain collimated, divergent or focused beam at the output.



Specifications

Lens material	AR coated Fused Silica Lenses
Screw Size	M22x0.75

Wavelength, nm	Expansion ratio	Beam expander size L, mm	Transmission, %	Catalogue number	Price, EUR
1064	2X	51	>96	160-0021	235
1064	2.5X	51	>96	160-0251	235
1064	3X	68	>96	160-0031	235
1064	4X	75	>96	160-0041	235
1064	5X	73	>96	160-0051	235
1064	6X	75	>96	160-0061	235
1064	8X	77	>96	160-0081	235
1064	10X	70	>96	160-0101	235
532	2X	51	>96	160-0022	235
532	2.5X	51	>96	160-0252	235
532	3X	68	>96	160-0032	235
532	4X	75	>96	160-0042	235
532	5X	73	>96	160-0052	235
532	6X	75	>96	160-0062	235
532	8X	77	>96	160-0082	235
532	10X	70	>96	160-0102	235
355	4X	75	>96	160-0043	250
355	6X	75	>96	160-0063	250
355	8X	68	>96	160-0083	250
355	10X	71	>96	160-0103	250

Compact beam expanders of other expansion ratio are available upon request.

Related Product

Large Rod Small Mounting
Clamp (aluminium) 810-0062A
Find more at EksmaOptics.com



ZOOM BEAM EXPANDER

Features

- Adjustable 1X – 8X or 2X – 8X expansion ratio
- Adjustable divergence
- Galilean design

Compact Galilean type zoom beam expanders are designed for Nd:YAG fundamental and harmonic wavelengths: 1064 nm, 532 nm and 355 nm. Zoom beam expanders provide

1X – 8X or 2X – 8X continuous magnification with adjustable focus to correct for laser beam divergence.

Wavelength, nm	Expansion ratio	Input Clear Aperture, mm	Output Clear Aperture, mm	Length, mm	Catalogue number	Price, EUR
1064	1x-8x	12	33	162	165-1181	860
1064	2x-8x	12	33	143.3	165-1281	860
532	1x-8x	12	33	162	165-1185	860
532	2x-8x	12	33	139.9	165-1285	860
355	1x-8x	12	33	162	165-1183	1120
355	2x-8x	12	33	158.5	165-1283	860



Visit our e-shop www.eksmaoptics.com and find the drawings of all zoom beam expanders.

Related Product

Universal Adjustable
Optics Mount 830-0035
Find more at EksmaOptics.com




TOP HAT BEAM SHAPING LENSES – FBS

Features

- Transforms Gaussian beams into flat-top beams
- Different shapes of Top Hat spots - round, square, line
- Easy to integrate into existing beam paths

FBS Series beam shapers are designed to transform collimated Gaussian beams into small and homogeneous Top Hat spots of square, round or line shapes. FBS beam shapers should be used in setup with focusing optics, and the working distance of a Top Hat beam shaper is determined by the focal length of the focusing optics.




Specifications

Material	UV Fused Silica
Transmission	>99%
Diameter	25.4 mm
Thickness	3 mm
Laser Damage Threshold	10 J/cm ² @ 1064 nm, 10 ns 5 J/cm ² @ 532 nm, 10 ns 3 J/cm ² @ 355 nm, 10 ns


Operation Requirements

Input beam	Gaussian beam (TEM ₀₀), M ² < 1.4
Input beam diameter	Fixed, ±5% tolerance
Operation wavelength	Fixed
Optical setup	Clear apertures along beam path at least 2.2x larger than the beam size @ 1/e ²


Without FBS Beam Shaper




Gaussian-profile at focal plane




Diffraction limited Gaussian profile




With FBS Beam Shaper




Top-Hat-profile at focal plane



Near diffraction limited Top Hat profile




FBS2 – Square Top Hat Profile



Spot Geometry

Top Hat width	approximately $2 \times \lambda \times (f/d)$, with f = focal length, d = beam diameter @ $1/e^2$
Efficiency	up to 90%
Homogeneity	ca. $\pm 2.5\%$ (rel. to average intensity of the Top Hat plateau)
Side modes (strongest)	$\sim 16.5\times$ weaker than line plateau ($< 1.5\%$ of input energy)
Depth of focus (DOF)	$\sim 60\%$ of the Rayleigh length


FBSR – Round Top Hat Profile



Spot Geometry

Top Hat diameter	approximately $2 \times \lambda \times (f/d)$, with f = focal length, d = beam diameter @ $1/e^2$
Efficiency	up to 95%
Homogeneity	ca. $\pm 2.5\%$ (rel. to average intensity of the Top Hat plateau)
Side modes (strongest)	$\sim 70\times$ weaker than line plateau ($< 1.5\%$ of input energy)
Depth of focus (DOF)	$\sim 30\%$ of the Rayleigh length

FBSL – Line Top Hat Profile



Spot Geometry

Line length (Top Hat)	approximately $2 \times \lambda \times (f/d)$, with f = focal length, d = beam diameter @ $1/e^2$
Line width (Gaussian)	Similar to the diameter of Gaussian spot in the same optical configuration
Efficiency	up to 92.5%
Homogeneity	ca. $\pm 2.5\%$ (rel. to average intensity of the Top Hat plateau)
Side modes (strongest)	$\sim 15\times$ weaker than line plateau ($< 1.7\%$ of input energy)
Depth of focus (DOF)	$\sim 50\%$ of the Rayleigh length

Ordering information

FBS2-1064-1.0

Model name:

- FBS2 – Square profile
- FBSR – Round profile
- FBSL – Line profile

Operation wavelength:


- 1064 nm
- 1030 nm
- 532 nm
- 515 nm
- 355 nm
- 343 nm

Input beam diameter (@ $1/e^2$):

- From 1.0 to 6.0 mm,
increment: 0.5 mm

There are different possibilities to integrate the FBS beam shaper into an optical setup.

1. Beam shaper directly in front of a focusing optic/objective




By introducing the FBS beam shaper into the beam path in front of a lens/objective the initial diffraction limited Gaussian spot will be transformed into a homogeneous Top-Hat profile.

When a Gaussian TEM₀₀ input beam with a diameter of 5 mm @ 1/e² is used the diameter of the free aperture along the total beam path have to be at least 11 mm (better 13 mm).

If for example a wavelength with 532 nm, a Gaussian TEM₀₀ input beam with a diameter of 5 mm @ 1/e² and a focusing lens with f=160 mm is used, ones will get a homogeneous Top Hat profile with a diameter of 34 µm.


2. Beam shaper in front of a beam expander



There is also the possibility to introduce the FBS beam shaper into the beam path in front of a beam expander. This leads to a higher numerical aperture of the focused beam and to a smaller Top Hat profile.


Example: A Gaussian beam with a diameter of 5 mm @ 1/e² illuminates the FBS beam shaper and is afterwards increased by a beam expander to a beam diameter of 8 mm. With an focusing optic with f=50 mm the user can generate a Top Hat with a diameter of 7 µm. The needed free aperture increases with the expanded beam. For a beam with a diameter of 8 mm the free aperture has to be at least 18 mm.

3. Beam shaper within a beam expander




A further and even more flexible possibility is to introduce the FBS beam shaper into the beam path within a beam expander. The user has the possibility for an easy "fine tuning" of beam diameter at the position of FBS beam shaper by shifting shaper along z-axis.

Scribing of CIGS-solar cells




- Wasted area, reducing efficiency → need of smallest scribing lines
- Cut quality influence efficiency → need of small HAZ, no debris, smooth edges
- High scanning speed for high throughput → need of small pulse overlap

P1 – „Scribing“




Gaussian Profile




FBS-Top-Hat Profile
small overlap, smooth edges

Removal of a front contact in ZnO(1 µm)/CIGS/Mo/PI structure. Laser PL10100/SH, 10 ps, 370 mW, 100 kHz, 532 nm; scanning speed 4.3 m/s, single pass.

P3 – „Scribing“



Gaussian Profile



FBS-Top-Hat Profile
small HAZ, smooth edges

Tilted SEM pictures of the P3 scribe in ZnO(1 µm)/CIGS/ Mo/PI structure. Laser PL10100/SH, 10 ps, 370 mW, 100 kHz, 532 nm; scanning speed 60 mm/s, single pass.

Raciukaitis et. al, JLMN-Vol. 6, No. 1, 2011

Recommended Accessories

Zoom Beam Expander
See page 5.4



Two Axes Translation Polarizer Holder
840-0240
Find more at
EksmaOptics.com



CONTINUOUSLY VARIABLE ATTENUATOR / BEAMSPLITTER – 990-0060

Features

- Divides laser beam into two beams of manually adjustable intensity ratio
- Convenient 90° angle between reflected and transmitted beams
- Negligible beam deviation
- Large dynamic range
- Broadband transmission
- Weight – 0.16 kg



Continuously Variable Attenuator/ Beamsplitter is designed to be used for laser pulses as short as 100 fs. It consists of 2 high-performance polarizing optics components placed in precision opto-mechanical holder 840-0197. Variable attenuator/beamsplitter incorporates a high-performance Polarizing Cube Beamsplitter which reflects s-polarized light at 90° while transmitting p-polarized light.

A rotating $\lambda/2$ waveplate is placed in the incident polarized laser beam. The intensity ratio of those two beams may be continuously varied without alteration of other beam parameters by rotating the waveplate. The intensity of either exit beam, and their intensity ratio, can be controlled over a wide dynamic range. Pure p-polarization could be selected for maximum transmission, or pure s-polarization for maximum attenuation of the transmitted beam.

ACHROMATIC AIR-SPACED WAVEPLATE AND HIGH POWER BROADBAND CUBE POLARIZING BEAMSPLITTER

Specifications

Extinction ratio	Ts/Tp < 1:200
Clear aperture	11 mm

for Broadband Region

Central wavelength, nm	LDT, J/cm ²	Catalogue number	Price, EUR
450-680	1 ¹⁾	990-0060-11VIS	1060
700-1000	2 ²⁾	990-0060-11IR	1060

¹⁾ LDT measured at 532 nm, 10 Hz, 10 ns pulses.

²⁾ LDT measured at 1064 nm, 10 Hz, 10 ns pulses.

MULTIPLE ORDER HALF WAVEPLATE AND HIGH POWER CUBE POLARIZING BEAMSPLITTER

Specifications

Extinction ratio	Ts/Tp < 1:500
Clear aperture	11 mm

Central wavelength, nm	LDT, J/cm ² *	Catalogue number	Price, EUR
1064	15	990-0061-11	740
1030	15	990-0062-11	740
800	8	990-0063-11	740
532	6	990-0064-11	740
355	3	990-0065-11	770

* LDT measured at designed wavelength, 10 Hz, 10 ns pulses.

VARIABLE ATTENUATORS FOR LINEARLY POLARIZED LASER BEAM – 990-0070

Features


- Divides laser beam into two parallel beams of manually adjustable intensity ratio
- Large dynamic range
- Transmitted beam shift ~ 0.5 mm
- High optical damage threshold




Note: Movable base **820-0090**, Rod Holder **820-0050-02** and standard rod should be ordered separately.

This variable attenuator/beamsplitter consists of special design opto-mechanical Adapter and precision opto-mechanical holder 840-0197. Two Thin Film Brewster type polarizers, which reflect s-polarized light while transmitting p-polarized light, are housed into Adapter. Quartz Half Waveplates are housed in rotating holder 840-0197.

The intensity ratio of those two beams may be continuously varied without alteration of other beam parameters by rotating the waveplate. The intensity of either exit beam,



or their intensity ratio, can be controlled over a wide dynamic range. P-polarization could be selected for maximum transmission, or high-purity s-polarization could be reflected when maximum attenuation of the transmitted beam takes place. The holder 840-0197 allows to adjust Angle Of Incidence of the Thin Film Brewster type polarizers by $\pm 2^\circ$ and to get the maximum polarization contrast.



For Nd:YAG Laser Applications

Aperture diameter	17 mm
Damage threshold	5 J/cm ² pulsed at 1064 nm, typical
Polarization Contrast (after 1st polarizer)	>1:200
Polarization Contrast (after 2nd polarizer)	>1:500
Weight	0.35 kg

For Nd:YAG Laser Applications

Wavelength, nm	Catalogue number	Price, EUR
266	990-0070-266H *	1050
355	990-0070-355	780
532	990-0070-532	680
1064	990-0070-1064	680

Multi order half waveplate is housed in rotating holder 840-0197 for Nd:YAG laser pulses (laser damage threshold: 5 J/cm² pulsed at 1064 nm, typical).

* With Zero Order Air-Spaced half waveplate.

For Femtosecond Applications

Aperture diameter	17 mm
Damage threshold	>10 mJ/cm ² , 50 fs pulse at 800 nm, typical
for high power laser applications	>100 mJ/cm ² , 50 fs pulse at 800 nm, typical
Time dispersion	t<4 fs for 100 fs Ti:Sapphire laser pulses
Polarization Contrast (after 1st polarizer)	>1:200
Polarization Contrast (after 2nd polarizer)	>1:500
Weight	0.35 kg

For Femtosecond Applications

Wavelength, nm	Catalogue number	Price, EUR
257	990-0070-257	975
266	990-0070-266	975
343	990-0070-343	870
400	990-0070-400	770
390-410	990-0070-400B	920
515	990-0070-515	770
505-525	990-0070-515B	920
800	990-0070-800	770
780-820	990-0070-800B	920
1030	990-0070-1030	770
1010-1050	990-0070-1030B	920

Zero order optically contacted half waveplate is housed in rotating holder 840-0197 for femtosecond laser pulses (laser damage threshold: >10 mJ/cm², 50 fsec pulse, 800 nm typical).

For High Power Femtosecond Laser Applications

Wavelength, nm	Catalogue number	Price, EUR
257	990-0070-257H	1050
266	990-0070-266H	1050
343	990-0070-343H	945
400	990-0070-400H	845
390-410	990-0070-400HB	995
515	990-0070-515H	845
505-525	990-0070-515HB	995
800	990-0070-800H	845
780-820	990-0070-800HB	995
1030	990-0070-1030H	845
1010-1050	990-0070-1030HB	995

Zero Order Air-Spaced half waveplate is housed in rotating holder 840-0197 for high power femtosecond applications (laser damage threshold: >100 mJ/cm², 50 fsec pulse, 800 nm typical).


Related Products

Beam dumps
990-0800,
990-0820

See page 5.19



MOTORIZED VARIABLE ATTENUATOR FOR LINEARLY POLARIZED LASER BEAM – 990-0070M



This motorized variable attenuator/beamsplitter consists of special design opto-mechanical Adapter and precision opto-mechanical holder 840-0193. Two Thin Film Brewster type polarizers, which reflect s-polarized light while transmitting p-polarized light, are housed into Adapter. Quartz Half Waveplates are housed in motorized rotation stage 960-0161.

The intensity ratio of those two beams may be continuously varied without alteration of other beam parameters by rotating the waveplate. The intensity of either exit beam, or their intensity ratio, can be controlled over a wide dynamic range. P-polarization could be selected for maximum transmission, or high-purity s-polarization could be reflected when maximum attenuation of the transmitted beam takes place. The holder 840-0193 allows to adjust Angle Of Incidence of the Thin Film Brewster type polarizers by $\pm 2^\circ$ and to get the maximum polarization contrast.

Ordering information

Please note: these motorized variable attenuators for linearly polarized laser beam are provided without controller and power supply. If you would like to order the complete solution (controller 980-1045 and power supply: PS12-1.5-4), please add CP to code and 695 EUR to price.


Example:

990-0070-266M – motorized attenuator without controller and power supply.

Price – 1830 EUR

990-0070-266M+CP – motorized attenuator with controller and power supply.

Price – 2525 EUR



For Nd:YAG Laser Applications

Wavelength, nm	Catalogue number	Price, EUR
266	990-0070-266HM *	1830
355	990-0070-355M	1560
532	990-0070-532M	1460
1064	990-0070-1064M	1460

Multi order half waveplate is housed in Motorized Rotation Stage 960-0161 and Polarizer with adapter in Kinematic Optical Mount 840-0193 for Nd:YAG laser application (laser damage threshold: 5 J/cm², 10 ns pulses, 10 Hz at 1064 nm, typical).

* With Zero Order Air-Spaced half waveplate.

For Femtosecond Applications

Wavelength, nm	Catalogue number	Price, EUR
257	990-0070-257M	1755
266	990-0070-266M	1755
343	990-0070-343M	1650
400	990-0070-400M	1550
390-410	990-0070-400BM	1700
515	990-0070-515M	1550
505-525	990-0070-515BM	1700
800	990-0070-800M	1550
780-820	990-0070-800BM	1700
1030	990-0070-1030M	1550
1010-1050	990-0070-1030BM	1700

Zero order optically contacted half waveplate is housed in Motorized Rotation Stage 960-0161 and Polarizer with adapter in Kinematic Optical Mount 840-0193 for femtosecond laser application (laser damage threshold: >10 mJ/cm², 50 fsec pulse, 800 nm typical).

For High Power Femtosecond Applications

Wavelength, nm	Catalogue number	Price, EUR
257	990-0070-257HM	1830
266	990-0070-266HM	1830
343	990-0070-343HM	1725
400	990-0070-400HM	1625
390-410	990-0070-400HBM	1775
515	990-0070-515HM	1625
505-525	990-0070-515HBM	1775
800	990-0070-800HM	1625
780-820	990-0070-800HBM	1775
1030	990-0070-1030HM	1625
1010-1050	990-0070-1030HBM	1775

Zero Order Air-Spaced half waveplate is housed in Motorized Rotation Stage 960-0161 and Polarizer with adapter in Kinematic Optical Mount 840-0193 for high power femtosecond laser application (laser damage threshold: >100 mJ/cm², 50 fsec pulse, 800 nm typical).

BROADBAND VARIABLE ATTENUATOR FOR FEMTOSECOND LASER PULSES – 990-0070HBBi70

Features

- Divides laser beam into two parallel beams of manually adjustable intensity ratio
- Large dynamic range
- Transmitted beam shift ~ 2.6 mm
- High optical damage threshold




990-0070-800HBBi70



990-0070-800HBBi70M

This variable attenuator/beamsplitter consists of a special design opto-mechanical adapter and a precision opto-mechanical holder 840-0197. Two thin film polarizers, operating at AOI=70° and reflecting s-polarized light while transmitting p-polarized light, are housed into the adapter. A quartz zero order air-spaced half waveplate is housed into the rotating holder 840-0197.


The intensity ratio of outgoing two parallel beams may be continuously varied without alteration of other beam parameters by rotating the waveplate. The intensity of the



exit beam or outgoing beams intensity ratio can be controlled over a wide dynamic range.

P-polarized beam is transmitted straightly with a 2.6 mm shift and s-polarized beam (after 2 reflections) is parallel to the outgoing p-polarized beam, just separated by 28 mm.

The 840-0197 holder allows to adjust angle of incidence of the thin film polarizers by $\pm 2^\circ$ and to achieve the maximum polarization contrast.



Specifications

Aperture diameter	12 mm
Operating bandwidth	100 nm
Damage threshold	50 mJ/cm ² pulsed at 800 nm, 50 fsec, 50 Hz
Polarization contrast (after 1st polarizer)	>1:200
Polarization contrast (after 2nd polarizer)	>1:500

Manual attenuators

Wavelength, nm	Catalogue number	Price, EUR
750-850	990-0070-800HBBi70	1300
980-1080	990-0070-1030HBBi70	1300

Motorized attenuators

Wavelength, nm	Catalogue number	Price, EUR
750-850	990-0070-800HBBi70M	2080
980-1080	990-0070-1030HBBi70M	2080

Ordering information

Please note: these motorized variable attenuators for linearly polarized laser beam are provided without controller and power supply. If you would like to order the complete solution (controller 980-1045 and power supply: PS12-1.5-4), please add CP to code and 695 EUR to price.



Example:

990-0070-800HBBi70 – motorized attenuator without controller and power supply.



Price – 2080 EUR

990-0070-800HBBi70+CP – motorized attenuator with controller and power supply.


Price – 2775 EUR



990-0070-800HBBi70



990-0070-800HBBi70M



VARIABLE ATTENUATORS FOR LINEARLY POLARIZED LASER BEAM – 990-0071

Features


- Divides laser beam into separated by 68° angle two beams of manually adjustable intensity ratio
- Large dynamic range
- Transmitted beam shift ~0.5 mm
- High Optical damage threshold




Note: Solid Base Height Extender **820-0210** and Standard Rod **820-0020-20** should be ordered separately

This variable attenuator/beamsplitter consists of special design opto-mechanical adapter for polarizer at 56° 840-0117A or 840-0118A and precision opto-mechanical holder 840-0197. Thin Film Brewster type polarizer, which reflect s-polarized light at 56° while transmitting p-polarized light, is housed into adapter for polarizer at 56°. Quartz Half Waveplates are housed in rotating holder 840-0197.

The intensity ratio of those two beams may be continuously varied without alteration of other beam parameters by rotating the waveplate. The intensity of either exit beam,



or their intensity ratio, can be controlled over a wide dynamic range. P-polarization could be selected for maximum transmission, or high-purity s-polarization could be reflected when maximum attenuation of the transmitted beam takes place. The holder 840-0197 allows to adjust Angle Of Incidence of the Thin Film Brewster type polarizer by $\pm 2^\circ$ and to get the maximum polarization contrast.



For Nd:YAG Laser Applications

Aperture diameter	10 mm
Damage threshold	5 J/cm ² pulsed at 1064 nm, typical
Polarization Contrast	>1:200
Weight	0.25 kg

For Nd:YAG Laser Applications

Wavelength, nm	Catalogue number	Price, EUR
266	990-0071-266H *	690
355	990-0071-355	475
532	990-0071-532	475
1064	990-0071-1064	475

Multi order half waveplate is housed in rotating holder 840-0197 for Nd:YAG laser pulses (laser damage threshold: 5 J/cm² pulsed at 1064 nm, typical).

* With Zero Order Air-Spaced half waveplate.

For Femtosecond Applications

Aperture diameter	10 mm
Damage threshold	>10 mJ/cm ² , 50 fs pulse at 800 nm, typical
for high power laser applications	>100 mJ/cm ² , 50 fsec pulse, 800 nm typical
Time dispersion	t<4 fs for 100 fs Ti:Sapphire laser pulses
Polarization Contrast	>1:200
Weight	0.25 kg

For Femtosecond Applications

Wavelength, nm	Catalogue number	Price, EUR
257	990-0071-257	655
266	990-0071-266	655
343	990-0071-343	630
400	990-0071-400	580
390-410	990-0071-400B	680
515	990-0071-515	580
505-525	990-0071-515B	680
800	990-0071-800	580
780-820	990-0071-800B	680
1030	990-0071-1030	580
1010-1050	990-0071-1030B	680


Zero order optically contacted half waveplate is housed in rotating holder 840-0197 for femtosecond laser pulses (laser damage threshold: >10 mJ/cm², 50 fs pulse at 800 nm, typical).

For High Power Femtosecond Laser Applications

Wavelength, nm	Catalogue number	Price, EUR
257	990-0071-257H	720
266	990-0071-266H	720
343	990-0071-343H	695
400	990-0071-400H	645
390-410	990-0071-400HB	745
515	990-0071-515H	645
505-525	990-0071-515HB	745
800	990-0071-800H	645
780-820	990-0071-800HB	745
1030	990-0071-1030H	645
1010-1050	990-0071-1030HB	745

Zero Order Air-Spaced half waveplate is housed in rotating holder 840-0197 for high power femtosecond applications (laser damage threshold: >100 mJ/cm², 50 fsec pulse, 800 nm typical).

MOTORIZED VARIABLE ATTENUATOR FOR LINEARLY POLARIZED LASER BEAM – 990-0071M



This motorized variable attenuator/beamsplitter consists of special design opto-mechanical adapter for polarizer at 56° 840-0117A or 840-0118A and precision opto-mechanical holder 840-0193. Thin Film Brewster type polarizer, which reflect s-polarized light at 56° while transmitting p-polarized light, is housed into adapter for polarizer at 56°. Quartz Half Waveplates are housed in motorized rotation stage 960-0161.

The intensity ratio of those two beams may be continuously varied without alteration of other beam parameters by rotating the waveplate. The intensity of either exit beam, or their intensity ratio, can be controlled over a wide dynamic range. P-polarization could be selected for maximum transmission, or high-purity s-polarization could be reflected when maximum attenuation of the transmitted beam takes place. The holder 840-0193 allows to adjust Angle of Incidence of the Thin Film Brewster type polarizer by $\pm 2^\circ$ and to get the maximum polarization contrast.

Ordering information

Please note: these motorized variable attenuators for linearly polarized laser beam are provided without controller and power supply. If you would like to order the complete solution (controller 980-1045 and power supply: PS12-1.5-4), please add CP to code and 695 EUR to price.


Example:

990-0071-266M – motorized attenuator without controller and power supply.

Price – 1435 EUR

990-0071-266M+CP – motorized attenuator with controller and power supply.

Price – 2130 EUR



For Nd:YAG Laser Applications

Wavelength, nm	Catalogue number	Price, EUR
266	990-0071-266HM *	1500
355	990-0071-355M	1290
532	990-0071-532M	1260
1064	990-0071-1064M	1260

Multi order half waveplate is housed in Motorized Rotation Stage 960-0161 and Polarizer with adapter in Kinematic Optical Mount 840-0193 for Nd:YAG laser application (laser damage threshold: 5 J/cm², 10 ns pulses, 10 Hz at 1064 nm, typical).

* With Zero Order Air-Spaced half waveplate.

For Femtosecond Applications

Wavelength, nm	Catalogue number	Price, EUR
266	990-0071-266M	1435
343	990-0071-343M	1410
400	990-0071-400M	1360
390 – 410	990-0071-400BM	1460
515	990-0071-515M	1360
505 – 525	990-0071-515BM	1460
800	990-0071-800M	1360
780 – 820	990-0071-800BM	1460
1030	990-0071-1030M	1360
1010 – 1050	990-0071-1030BM	1460

Zero order optically contacted half waveplate is housed in Motorized Rotation Stage 960-0161 and Polarizer with adapter in Kinematic Optical Mount 840-0193 for femtosecond laser application (laser damage threshold: >10 mJ/cm², 50 fsec pulse, 800 nm typical).

For High Power Femtosecond Laser Applications

Wavelength, nm	Catalogue number	Price, EUR
266	990-0071-266HM	1500
343	990-0071-343HM	1475
400	990-0071-400HM	1425
390 – 410	990-0071-400HBM	1525
515	990-0071-515HM	1425
505 – 525	990-0071-515HBM	1525
800	990-0071-800HM	1425
780 – 820	990-0071-800HBM	1525
1030	990-0071-1030HM	1425
1010 – 1050	990-0071-1030HBM	1525

Zero Order Air-Spaced half waveplate is housed in Motorized Rotation Stage 960-0161 and Polarizer with adapter in Kinematic Optical Mount 840-0193 for high power femtosecond laser application (laser damage threshold: >100 mJ/cm², 50 fsec pulse, 800 nm typical).


VARIABLE ATTENUATOR FOR FEMTOSECOND LASER PULSES – 990-0072

Features

- Divides laser beam into two beams of manually adjustable intensity ratio separated by 68° angle
- Large dynamic range
- Transmitted beam shift ~1 mm
- High optical damage threshold



Check www.eksmaoptics.com for motorized version 990-0072M




This variable attenuator/beamsplitter consists of Polarizer Holder 840-0190-01 and Kinematic Mirror/Beamsplitter Mount 840-0056-12.

UVFS Thin Film Brewster type polarizer diameter 50.8 mm, which reflect s-polarized light while transmitting p-polarized light, is housed into Beamsplitter Mount 840-0056-12. A quartz Zero Order (optically contacted) Half Waveplate Ø25.4 mm (for femtosecond applications), quartz Zero Order Air-Spaced Half Waveplate (for high power femtosecond applications) or quartz Multi Order Half Waveplate Ø25.4 mm (for Nd:YAG laser applications) is housed in rotating polarizer holder 840-0180-A1 and placed in the incident linearly polarized laser beam.

The intensity ratio of those two separated and different polarized beams may be continuously varied without alteration of other beam parameters by rotating the waveplate. The intensity of either exit beam, or their intensity ratio, can be controlled over a wide dynamic range. P-polarization could be selected for maximum transmission, or high-

purity s-polarization could be reflected when maximum attenuation of the transmitted beam takes place.

The holder 840-0056-12 allows to adjust Angle Of Incidence of the Thin Film Brewster type polarizers by $\pm 4.5^\circ$ and to get the maximum extinction contrast. The mounts are on rods, rod holders and Movable Base 820-0090. The optical axis height from the table top can be adjusted in the range 78-88 mm. Other height can be offered as custom changing the standard rods and rod holders into higher.



For Nd:YAG Laser Applications

Clear Aperture diameter	22 mm
Damage threshold	>5 J/cm ² , 10 ns pulse, 10 Hz at 1064 nm, typical
Polarization Contrast	>1:200
Transmitted beam shift	~1 mm
Weight	0.45 kg

A quartz Multi Order Half Waveplate Ø25.4 mm is housed in rotating holder 840-0180-A1.

For Nd:YAG Laser Applications

Wavelength, nm	Catalogue number	Price, EUR
266	990-0072-266H *	1115
355	990-0072-355	795
532	990-0072-532	765
1064	990-0072-1064	785

* A quartz Zero Order Air-Spaced Half Waveplate clear aperture Ø22 mm is housed in rotating holder 840-0190-01.

For Femtosecond Applications

Clear Aperture diameter	22 mm
Damage threshold	>10 mJ/cm ² , 50 fs pulse at 800 nm, typical
for high power applications	>100 mJ/cm ² , 50 fs pulse at 800 nm, typical
Polarization Contrast	>1:200
Transmitted beam shift	~1 mm
Weight	0.45 kg

A quartz Zero Order (optically contacted) Half Waveplate (for femtosecond applications) or Zero Order Air-Spaced Half Waveplate (for high power applications) Ø25.4 mm are housed in rotating holder 840-0190-01.

For Femtosecond Applications

Wavelength, nm	Catalogue number	Price, EUR
266	990-0072-266	980
343	990-0072-343	925
400	990-0072-400	895
515	990-0072-515	895
800	990-0072-800	910
780 – 820	990-0072-800B	1010
1030	990-0072-1030	920
1010 – 1050	990-0072-1030B	1010

For High Power Femtosecond Laser Applications

Wavelength, nm	Catalogue number	Price, EUR
266	990-0072-266H	1115
343	990-0072-343H	1060
400	990-0072-400H	1030
515	990-0072-515H	1030
800	990-0072-800H	1045
780 – 820	990-0072-800HB	1145
1030	990-0072-1030H	1055
1010 – 1050	990-0072-1030HB	1145

VARIABLE ATTENUATOR FOR FEMTOSECOND AND Nd:YAG LASER PULSES – 990-0073



Features

- Divides laser beam into two beams of manually adjustable intensity ratio separated by 68° angle
- Large dynamic range
- Transmitted beam shift ~1.4 mm
- High optical damage threshold
- Motorized version available on request




This variable attenuator/beamsplitter consists of Polarizer Holder 840-0180-A2 and Kinematic Mirror/Beamsplitter Mount 840-0056-13. UVFS Thin Film Brewster type polarizer Ø76.2 mm, which reflect s-polarized light while transmitting p-polarized light, is housed into Beamsplitter Mount 840-0056-13. A quartz Zero Order (optically contacted) Half Waveplate Ø40 mm (for femtosecond applications), Zero Order Air-Spaced Half Waveplate (for high power femtosecond applications) or quartz Multi Order Half Waveplate Ø40 mm (for Nd:YAG laser applications) is housed in rotating polarizer holder 840-0180-A2 and placed in the incident linearly polarized laser beam.

The intensity ratio of those two separated and different polarized beams may be continuously varied without alteration of other beam parameters by rotating the waveplate. The intensity of either exit beam,



or their intensity ratio, can be controlled over a wide dynamic range. P-polarization could be selected for maximum transmission, or high-purity s-polarization could be reflected when maximum attenuation of the transmitted beam takes place.

The holder 840-0056-13 allows to adjust Angle Of Incidence of the Thin Film Brewster type polarizers by $\pm 4.5^\circ$ and to get the maximum extinction contrast. The mounts are on rods, rod holders and Movable Base 820-0090. The optical axis height from the table top can be adjusted in the range 92-98 mm. Other height can be offered as custom changing the standard rods and rod holders into higher.



For Nd:YAG Laser Applications

Clear Aperture diameter	36 mm
Damage threshold	>5 J/cm ² , 10 ns pulse, 10 Hz at 1064 nm, typical
Polarization Contrast	>1:200
Transmitted beam shift	~ 1.4 mm
Weight	0.6 kg

Quartz Multi Order Half Waveplate Ø40 mm is housed in rotating polarizer holder 840-0180-A2.

For Nd:YAG Laser Applications

Wavelength, nm	Catalogue number	Price, EUR
266	990-0073-266H *	1820
355	990-0073-355	1490
532	990-0073-532	1470
1064	990-0073-1064	1545

* Zero Order Air-Spaced half waveplate is housed in rotating holder.

For Femtosecond Applications

Clear Aperture diameter	36 mm
Damage threshold	>10 mJ/cm ² , 50 fs pulse at 800 nm, typical
for high power applications	>100 mJ/cm ² , 50 fs pulse at 800 nm, typical
Polarization Contrast	>1:200
Transmitted beam shift	~ 1.4 mm
Weight	0.6 kg

A quartz Zero Order (optically contacted) Half Waveplate Ø40 mm (for femtosecond applications) or Zero Order Air-Spaced Half Waveplate (for high power applications) is housed in rotating polarizer holder 840-0180-A2.

For Femtosecond Applications

Wavelength, nm	Catalogue number	Price, EUR
266	990-0073-266	1720
343	990-0073-343	1590
400	990-0073-400	1570
515	990-0073-515	1570
800	990-0073-800	1590
780 – 820	990-0073-800B	1820
1030	990-0073-1030	1645
1010 – 1050	990-0073-1030B	1880

For High Power Femtosecond Laser Applications

Wavelength, nm	Catalogue number	Price, EUR
266	990-0073-266H	1820
343	990-0073-343H	1690
400	990-0073-400H	1670
515	990-0073-515H	1670
800	990-0073-800H	1690
780 – 820	990-0073-800HB	1920
1030	990-0073-1030H	1745
1010 – 1050	990-0073-1030HB	1980

COMPACT MOTORIZED LASER POWER ATTENUATOR – 990-0075



990-0075 series laser power attenuator is a compact motorized device for laser power control of linearly polarized beam. Its operation principle consists of changing laser beam polarization by rotation of zero-order half-waveplate and then passing the beam through a fixed Brewster type thin film polarizer which transmits p-polarized and reflects s-polarized laser beam. The intensity ratio of s- and p- polarized beams can be continuously varied without alteration of other beam parameters by motorized rotation of the waveplate. The attenuation level is controlled by the software in 0.1 – 97% range.

The device combines unique mechanical design which ensures repeatability and high stability of performance. All optical components of the attenuator have high LIDT coatings and provide stable and reliable performance even using device with high power lasers in industrial applications.

A secondary laser beam with s-polarization can be rejected from the laser power attenuator unit to an external beam dump or utilized for the particular application. A standard compact external beam dump is optionally offered with this attenuator and is suitable for lasers with average power up to 6 W. This beam dump stops secondary s-polarized beam in the attenuator and allows to avoid any thermal effects or stress in the housing of the device.

Features

- Compact design
- High optical damage threshold
- Full solution - includes controller, software, power supply and USB cable
- Standard models for the most popular laser wavelengths are offered ex-stock

Standard Kit includes:

- Motorized laser power attenuator
- Controller
- Software
- Power supply (DC 12 V)
- USB cable (1.5 m)

Optical Specifications


Clear input aperture	Ø18 mm
Clear output aperture	Ø18 mm
Power attenuation range	<0.1 % to >98 %
LIDT coating	>10 J/cm ² , 10 ns @ 1064 nm
Close to open time	< 0.2 sec
Resolution	175,542 µsteps in full rotation (0.002 deg, 7.2 arcsec, 0.035 mrad)
Accuracy (repeatability after 10,000 positions without homing)	±10 µsteps (±0.02 deg, less than ±0.035 %)
Motor	2 phase stepper motor, 200 steps with 256 µstepping

Mechanical Specifications

	Length	Width	Height
Attenuator	76 mm	36 mm	58 mm
Attenuator with beam dump	76 mm	52 mm	58 mm
Controller	125 mm	53 mm	31 mm

Operating Conditions

Operating temperature	10 to 40 °C
Storage temperature	-15 to 50 °C



Laser Power Attenuators

Wavelength, nm	LIDT	Catalogue number	Price, EUR
343	3 J/cm ² , 10 ns, 50 Hz @ 343 nm	990-0075-343M	1420
355	3 J/cm ² , 10 ns, 50 Hz @ 355 nm	990-0075-355M	1420
390 – 410	3 J/cm ² , 10 ns, 50 Hz @ 400 nm	990-0075-400M	1420
510 – 520	5 J/cm ² , 10 ns, 50 Hz @ 515 nm	990-0075-515M	1310
532	5 J/cm ² , 10 ns, 50 Hz @ 532 nm	990-0075-532M	1310
780 – 820	8 J/cm ² , 10 ns, 50 Hz @ 800 nm	990-0075-800M	1310
1020 – 1040	10 J/cm ² , 10 ns, 50 Hz @ 1030 nm	990-0075-1030M	1310
1064	10 J/cm ² , 10 ns, 50 Hz @ 1064 nm	990-0075-1064M	1310

Additional Accessories

Catalogue number	Description	Price, EUR
990-0075ABD	Attachable Beam Dump <6W with coated protective window	165
990-0075SBD	Separated Beam Dump < 30 W	175
990-0075C1	RS232 Cable, 1.8 m	15
990-0075C5	RS232 Cable, 5 m	20
990-0075C10	RS232 Cable, 10 m	30
990-0075C15	RS232 Cable, 15 m	35
990-0075RA	DIN35 Rail Adapter	15

COMPACT VARIABLE LASER POWER ATTENUATOR – 990-0076

Features

- Attenuation range for reflected beam: 2% – 99.9%
- Attenuation range for transmitted beam: 0.1% – 98%
- Convenient 90° angle between reflected and transmitted beams
- No beam shift at exit




990-0076 Variable Attenuator incorporates a high-performance polarizing cube beamsplitter which reflects s-polarized light at 90° while transmitting p-polarized light. A rotatable $\lambda/2$ waveplate is placed in the incident polarized laser beam.

The intensity ratio of the two outgoing beams may be continuously varied without alteration of other beam parameters by rotating the half waveplate. The intensity of either exit beam, or their intensity ratio, can be controlled over a wide dynamic range. Pure p-polarization could be selected for maximum transmission, or pure s-polarization for maximum reflection.

Specifications

Clear aperture diameter	10 mm
Polarization contrast	>1:1000
Weight	0.1 kg

Wavelength, nm	LIDT	Catalogue number	Price, EUR
343	3 J/cm ² @ 355 nm, 10 Hz, 10 ns	990-0076-343	825
355	3 J/cm ² @ 355 nm, 10 Hz, 10 ns	990-0076-355	825
515	6 J/cm ² @ 532 nm, 10 Hz, 10 ns	990-0076-515	800
532	6 J/cm ² @ 532 nm, 10 Hz, 10 ns	990-0076-532	800
1030	15 J/cm ² @ 1030 nm, 10 Hz, 10 ns	990-0076-1030	800
1064	15 J/cm ² @ 1064 nm, 10 Hz, 10 ns	990-0076-1064	800



FILTERS HOLDER WITH 90° FLIP – 990-0400

Features

- Allows stacking of 5 filters of Ø25.4 mm (1"), or 3 filters of Ø50.8 (2")
- Fast flipping in and out of beam path
- Available to be used in 90° position
- Has one M4, two M6 and two holes Ø6.4mm for mounting on posts or table bases
- Large aperture allows to attenuate large diameter laser beam
- Black Anodized Aluminium and Brass screws



990-0415



990-0423

The holder of 1 inch filters 990-0415 allows the fixation of up to 5 filters into 1 inch optics ring holders. The thickness of optical filters (or any other optical elements) to be held is from 0.5 mm to 8.0 mm. Filters can be easily replaced in ring holders. This filter holder allows fast filter removal from beam path flipping it at 90° position. Any position of filters can be fixed with fixing screw. The firm 0° position can be fixed with the second brass screw (included).

The holder of 2 inch filters 990-0423 allows the fixation of up to 3 filters into 2 inch optics ring holders. The thickness of optical filters (or any other optical elements) to be held is from 0.5 mm to 14.0 mm.

The holder 990-0415ND is the same holder 990-0415 but with Neutral Density filters that operates as step energy attenuator and allows adjusting transmission from 100% (all 5 filters are at 90° position) till 0.015% (all 5 filters are at 0° position) at visible region. If you need other adjustment you can choose any other Neutral Density filter Ø25.4 mm.

Using the holder 990-0415 with various color glass or dielectric filters various transmitted band pass regions can be achieved. The Filters Holder with 90° Flip is made of black anodized aluminium and brass screws.





Acceptable filters number	Suitable filters diameter, mm	Clear aperture diameter, mm	Weight, kg	Catalogue number	Price, EUR
5	25.4	23	0.16	990-0415	170
5	25.4	23	0.19	990-0415ND	265
3	50.8	48	0.22	990-0423	160



990-0415 at 0° position
(Note: Solid base height extender 820-0210 should be ordered separately)



990-0423 at 0° position
(Note: Solid base height extender should be ordered separately)



990-0415 at 0° or 90° position
(Note: Solid base height extender 820-0210 should be ordered separately)

Related Products

Neutral Density Filters Ø25.4 mm

See page 1.14

AIR-COOLED BEAM DUMP – 990-0800



990-0800

Beam Dump 990-0800 is designed to block CW or pulsed laser beams. It can be used on beams of up to 50 W in the wavelength range from 0.1 to 30 μm .


Due to the design of the beam dump, even if the non-reflective coating is damaged by high intensity pulses, there is no backward reflection.

990-0801 is a smaller beam dump designed to block a CW or a pulsed laser beam. It can be used on beams of up to 5 Watts in the wavelength range from 0.1 to 30 μm .


Specifications

Wavelength range	0.1 – 30 μm
Laser type	pulsed, CW

Aperture	Description	Catalogue number	Price, EUR
48 mm	for beams up to 50 W	990-0800	169
20 mm	for beams up to 5 W	990-0801	119



Drawing of 990-0800



Drawing of 990-0801

WATER-COOLED BEAM DUMP – 990-0820



990-0820

Beam Dump 990-0820 is designed to block CW or pulsed laser beams. It is mainly intended for beams 2 inch wide.

The dump is best suited for beams of up to 1 kW from 0.1 – 30 μm wavelength range. Even if the non-reflective coating is damaged by high intensity pulses, the beam is not reflected back into your optical scheme.

The dump mounts on M6 hole on its back.

Specifications

Wavelength range	0.1 – 30 μm
Max. handling power	1 kW
Max. energy	50 J (20 Hz)
Acceptance aperture	48 mm (1.89")
Laser type	pulsed, CW
Weight	1.2 kg

Catalogue number	Price, EUR
990-0820	289

